



## PRODUCT: DR200

For prolonged Holter monitoring periods, eCardio offers the DR200 Holter Monitor. The DR200's all digital recording eliminates mechanical malfunction and can store up to 14-days of data on a single flash card. The DR200 offers:

- All-digital recording eliminating mechanical malfunctions and increased reliability
- Removable digital memory (compact flash card) permitting immediate download and multiple use
- Ergonomic design: An easy-to-read LCD screen which displays actual patient ECG on all channels to ensure quality hook-up and lead placement.
- Multi-language capabilities: Recorder operation can be displayed in numerous languages.
- Pacemaker detection standard
- 3-Channel continuous monitoring for both pediatric and adult patients

### HOLTER SPECIFICATIONS

- Recording duration: 14 days with new alkaline batteries
- Recording bandwidth: 0.05 – 70 Hz
- Resolution: 12-bit, 12.50mV/LSB
- Sampling Rate: 180 samples per second
- Dimensions: 8.6 cm x 6.0 cm x 2.0 cm (3.9 in. x 2.4 in. x .8 in.)
- Weight: 99.3 grams (3.5 oz) including batteries and memory card
- Power: One (1) AA battery: alkaline, NiMH rechargeable or lithium
- Storage Capacity: 32 MB – 1 GB SD flashcard

