



## PRODUCT: POST-EVENT RECORDER

For transient symptomatic event detection, eCardio offers the post-symptomatic event recorder. The post-symptomatic event recorder is a lightweight, transtelephonic ECG recorder designed to document post-symptomatic cardiac activity. The monitor is designed for easy patient operation. The device offers:

- Simple operation with a single button for the manual recording of symptomatic events
- Lightweight, ergonomic features (2.5 oz) designed to improved ECG signal quality and enhance patient compliance
- Programmable memory nodes, advanced baseline correction and an auto-erase option
- Extended battery life, with more than 500 recordings and transmissions on a single battery, enhancing patient compliance

### INTENTIONS FOR USE

- Transient Symptomatic Event Detection
- Transient Ischemic Episodes
- ST Segment Analysis
- Detecting and Characterizing Symptomatic Arrhythmias
- Regulation of Anti-Arrhythmia Drug Therapy
- Drug Monitoring
- Post-myocardial infarction follow-up
- Cardiac Rehabilitation
- Variant Angina
- Pediatric Cardiology

